


STEALTH BELT®

OSTOMY SUPPORT BELTS

Double Ostomy Custom Worksheet



INSTRUCTIONS:

1. Take measurement A around the body above the highest stoma. Record measurement below.
2. Take measurement B around the body below the lowest stoma. Record the measurement below.
3. Take measurement C from the center of one stoma to the center of the other stoma. Record measurement below.
4. To find measurement D, find the center point of the highest stoma and follow the horizontal line down toward the other stoma. Once directly above center of the lowest stoma, measure down to the center point of that stoma. Record measurement below.
5. Indicate if you would like your belt in the horizontal style or vertical style for both bags? **LEFT-** _____ **RIGHT-** _____

LEFT Flange: _____ **Highest Stoma:** _____ **RIGHT Flange:** _____

A: _____

B: _____

C: _____

D: _____